



the art and science of smart patch testing™

## PG330: Tetrahydrofurfuryl methacrylat

**CAS#: 2455-24-5**

### Patient Information

Your patch test result indicates that you have a contact allergy to tetrahydrofurfuryl methacrylat. This contact allergy may cause your skin to react when it is exposed to this substance although it may take several days for the symptoms to appear. Typical symptoms include redness, swelling, itching, and fluid-filled blisters.

### Where is tetrahydrofurfuryl methacrylat found?

This substance is used to promote bonding in products. It is found in adhesives.

Avoid products that list any of the following names in the ingredients:

- 2-Propenoic acid, 2-methyl-, (tetrahydro-2-furanyl)methyl ester
- AI3-08497
- EINECS 219-529-5
- HSDB 5461
- Methacrylic acid tetrahydrofurfuryl ester
- NSC 32634
- SR 203
- Sartomer SR 203
- Tetrahydrofurfuryl methacrylate
- Tetrahydrofurfurylmethacrylate

### What are some products that may contain tetrahydrofurfuryl methacrylat?

#### Household Products

- Depend II Industrial Strength Adhesive
- Trim/Detail Adhesive (Kit with Activator)

\*For additional information about products that might contain **tetrahydrofurfuryl methacrylat**, go to the Household Product Database online (<http://householdproducts.nlm.nih.gov>) at the United States National Library of Medicine. These lists are brief and provide just a few examples. They are not comprehensive. Product formulations also change frequently. Read product labels carefully and talk to your doctor if you have any questions. These are general guidelines. Talk to your doctor for more specific instructions.